




Mechanical Ventilator Simulator



#### Compatibility

The first Ventilator paired with any patient simulators.



#### User-Friendly Interface

Simple, fast and effective user interaction with the simulator.



#### Mobility

View the instructor interface on any device with an internet browser.

**LungSim™** is a unique and immersive mechanical ventilator simulator that is able to be interfaced **with your human patient simulator**. This **model driven** software allows the user to be self-trained on the respiratory mechanisms (**standalone mode**) as well as to create advanced simulation scenarios on different patients with pulmonary diseases or acute respiratory failures when wirelessly linked to a patient simulator (**simulation mode**).

**LungSim™** incorporates real-time changes in vital signs, ventilatory parameters and more. Users can be trained

on define ventilator data, formulate diagnoses and make changes to best treat pathological patients. It allows instructors to simulate several patient conditions as for example chronic obstructive pulmonary diseases or acute respiratory distress syndromes simply changing patient parameters and ventilator values.

This software is composed by a realistic “user interface” with integrated controls that simulate a real ventilator and a friendly browser based “instructor interface” that can be used on PC, tablet, smartphone connected via lo-

cal network to the user interface and allows to use **LungSim™** during advanced simulation scenarios. The user interface incorporates built-in validated layouts in order to simulate different ventilators and situations. It allows the user to choose among 4 different ventilation modes:

- **Volume Control:** through this ventilation mode the user can define the tidal volume administered to the patient. The airway pressure results from the compliance of the lungs, airways resistances and the inhaled volume.
- **Pressure Control:** using this ventilation mode the user can define the inspiratory pressure which is administered to the patient's lungs.
- **Pressure Support:** using this assisted ventilation mode, that is patient-triggered (by pressure or airflow), the user can set the ventilator simulator to provide assistance when the patient makes a breathing effort. In this ventilation mode the trigger threshold represents a critical setting; actually an incorrect value of this parameter can produce the so called missed triggering phenomenon.
- **APRV:** this ventilation mode is a form of CPAP which utilises releases of the high pressure level to a low pressure level intermittently and so it allows the administration of inverse ratio ventilation. In this modality, the patient has the freedom to make unassisted spontaneous breaths.

The user interface allows also the user to directly change the ventilation settings (as the Positive End-Expiratory Pressure, the Respiratory Rate or the Tidal Volume) and to see on the virtual ventilator monitor several output parameters updating in real time (like Fraction of Inspired Oxygen, Insp : Esp Ratio, Expiratory Minute Volume, Inspiratory Minute Volume, Drive Pressure Value, etc.) and waveforms (like Pressure, Volume, Flow, Volume - Pressure loop, Flow - Pressure loop, Flow - Volume loop etc.).

The instructor is able to change the value of the parameters like the Lung Compliance, Thorax Compliance, Resistances, Respiratory Muscle Effort that describe the patient characteristics and this will clearly affects the Respiratory Rate, the Inspiration Time, the Oxygen saturation, the End-Tidal CO<sub>2</sub> and the Arterial Pressure of the

high-fidelity simulator.

Through the instructor interface, the instructor also have the possibility to load pathology presets – like the different steps of a COVID-19 typical case – or to simulate common adverse events that can typically occur during a real ventilation scenario – presence of water in the endotracheal tube, leaks in the tube or in the patient's lungs, presence of a mucus plug, or the presence of a bend in the tube.

**LungSim™ allows the user to train on the use of a mechanical ventilator without the need to have or buy a real mechanical ventilator.**

#### LungSim™ includes:

- LungSim case including a USB stick with LungSim software, lifetime license key and user manual
- Installation support and remote training
- Lifetime free software updates and upgrades

*Minimum System Requirements for the PC:*


*Display: 15" - Display Resolution: 1920 x 1080*

*CPU: Intel Core i-5 - RAM: 16 GB*

*Operating system: Windows 10 Pro (64 bit)*

*Language: English - SSD (Solid State Drive): 128 GB*

*Dedicated graphic card RAM: 2 GB*



*LungSim™ is a product from Accurate, a company blending international experience, scientific research, engineering and development of truly effective hi-tech educational solutions in the medical field.*

*Clinical use disclaimer: LungSim™ is licensed for use for educational purposes only. LungSim™ is not intended for clinical use. Accurate srl patent.*

FOR FURTHER INFORMATION [HTTP://LUNGSIM.ACCURATESOLUTIONS.EU/](http://LUNGSIM.ACCURATESOLUTIONS.EU/)